

Xylazine Drug Test Strip

PROCEDURE CARD

1

Test Liquid/Powder

- Set up the test in room temperature before use.
- Dilute the drug to be tested in water. The Center for Forensic Science Research and Education (CFSRE) recommends adding one scoop (5-10mg) of drug sample to 5mL of water. Refer to the advice of your local health or harm reduction authority on how much water and drug sample you should use.



Test Urine

- Urinate directly into the mouth of the urine container.




2

- Remove the test strip from its sealed pouch and use it as soon as possible. For best results, the test should be performed within one hour.



3

- Hold the strip by the end, where the product name (XYL) is printed. To avoid contamination, do not touch the strip membrane (the white section of the strip).
- Holding the strip vertically, dip the test strip and immerse the absorbent end in the liquid about 1/2 of the maximum mark line for at least 15 seconds.



Note: Do not immerse past the MAX line on the test strip.

4

- Remove the strip from the sample and place it on a non-absorbent flat surface.
- Start the timer and wait for the colored band(s) to appear. Read the test result at 5 minutes. **DO NOT READ RESULTS AFTER 10 MINUTES**

